

Sterile Purified Water, USP Grade

Product Number: WPW-USP-1X6

Lot Number: 20106162

Date of Manufacture: 06/2016

Expiration Date: 06/2019

Storage Temperature: 2 - 30° C

BASIS FOR RELEASE

CHEMICAL ANALYSIS			
Test	Method	Specification	Result
Appearance	Visual	Clear Liquid	Clear Liquid
Conductivity	USP <645>	<5.0 μS/cm	0.64 μS/cm
Sterility	USP <71>	No Growth	No Growth
TOC (before packaging)	USP <643>	<500 ppb	26.03 ppb

STATEMENT FOR USE

This product is for further manufacturing, laboratory or research use only.
Not intended for direct human or veterinary parenteral administration or where prohibited by law.

Approval

Quality Assurance- Print Name	Quality Assurance- Sign Name:	Date:
Susan Riggs	<i>Susan Riggs</i>	17 Jun 2016