

Bovine Serum Albumin, Fatty Acid Poor – 10 kg packages

Product number: RMB-002-203 Lot Number: 20120814A
 Date of Manufacture: 08/2012 Expiration date: 08/2017
 Storage temperature: 2 – 8°C

Not for human or animal use.

Basis for Release

Test / Method <i>Tests performed in 4% solution unless otherwise noted</i>	Result
Appearance	Light powder
pH (1% solution)	7.4
Purity	98.7%
Calcium	< 3.0 mg/dL
Chloride	< 10 mEq/L
Potassium	< 1.0 mEq/L
Sodium	10.0 mEq/L
Fatty Acid	<0.5 mg/g
Microbial	Passed
Moisture	2.2 %
Filtration	Solution sterile filtered through 0.2µm prior to lyophilization

Approval		
Quality Assurance – Print Name: <i>Christy Sims</i>	Quality Assurance – Sign Name: <i>CSims</i>	Date: <i>14 Sept 12</i>

Fatty Acid Poor Bovine Serum Albumin is derived from bovine serum, processed from blood collected at a USDA licensed establishment and donor animals inspected ante and post mortem at the abattoir by a U.S. veterinary services inspector.