

Bovine LDL Concentrate –SAMPLE

Product number: RMB-001-007 Lot Number: 20130307DL

Date of Manufacture: 03/2013

Storage temperature: 2 – 8°C

Not for human or animal consumption

Basis for Release

Test	Result
Cholesterol Concentration	2412 mg/dL
LDL	1014 mg/dL
HDL	860 mg/dL
Triglycerides	456 mg/dL
Sodium	<80 mEq/L
Potassium	6.0 mEq/L
Chloride	<60 mEq/L
Total Protein	9.2 g/dL
Iron	236 mEq/dL
Appearance	Opaque solution
Filtration	0.2 micron

Approval		
Quality Assurance – Print Name: Abbigail Miller	Quality Assurance – Sign Name: <i>Abbigail Miller</i>	Date: 03/22/2013

Bovine LDL Concentrate is derived from bovine serum, processed from blood collected from healthy animals, sourced from the US.