

**Bovine Serum Albumin, lyophilized powder**

Product number: RMB-002-101 Lot Number: 20130318A

Date of Manufacture: 03/2013 Expiration date: 03/2018

Storage temperature: 2 – 8°C

*For evaluation purposes only – not for human or animal consumption*

**Basis for Release**

Test / Method <i>Tests performed in 4% solution unless otherwise noted</i>	Result
pH	7.2
Solubility	pass
Protein, Total	3.9 g/dL
Purity	98%
Calcium	3.3 mg/dL
Chloride	< 25 mEq/L
Potassium	0.3 mEq/L
Sodium	< 20 mEq/L
Microbial	0 cfu/mL

Approval		
Quality Assurance – Print Name:	Quality Assurance – Sign Name:	Date:
<i>William D. Pfitter</i>	<i>William D. Pfitter</i>	<i>04/04/2013</i>

Bovine albumin is derived from bovine serum; processed from blood collected from healthy animals, sourced from GBR-1 rated countries. Characterization assays are performed on various lots of this product, and will be provided upon request.