

Bovine Serum Albumin-IgG Free

Product number: BSA-AAF-XXX Lot Number: 20141006AF
 Date of Manufacture: 10/2014 Expiration date: 10/2019
 Storage temperature: 2 – 8°C



Not for human or animal consumption.

Basis for Release

Test / Method <i>Tests performed in 4% solution unless otherwise noted</i>	Specification	Result
Origin	New Zealand	New Zealand
Appearance	Off White to Light Brown Crystal Powder	Off White to Light Brown Crystal Powder
pH (1% solution)	Check and Record	7.66
Purity	≥ 96.0%	100%
Calcium	≤ 7.5 mg/dL	< 5 mg/dL
Potassium	< 1.0 mEq/L	< 1 mEq/L
Sodium	< 50 mEq/L	< 5 mEq/L
Microbial	< 100 CFU/g	0 CFU/g
IgG content	< 10 ppm	0.8 ppm
Endotoxin	≤ 10 EU/mg	< 5.6 EU/mg

Approval

Quality Assurance – Print Name: <i>Stacey Black</i>	Quality Assurance – Sign Name: <i>Stacey Black</i>	Date: <i>08/26/2015</i>
--	---	----------------------------

Bovine albumin is derived from raw bovine serum processed from blood collected from healthy animals and sourced from Category A rated countries.