

Bovine Serum Albumin – 0.2 µm filtered, Polymer enhanced 30% Solution with 0.1% Sodium Azide

Product number: BSA-BSL-XXX Lot Number: 20160412LB

Date of Manufacture: 04/2016 Expiration date: 04/2019

Storage temperature: 2 – 8°C



Not for human or animal consumption

Basis for Release

Test/Method	Specification	Result
Origin	United States	Pass
Appearance	Light Brown to Amber Liquid	Pass
pH	Check and Record	7.2
Total Protein	As Reported	31 g/dL
Sodium	As Reported	95 mEq/L
Potassium	As Reported	1.0 mEq/L
Chloride	As Reported	40 mEq/L
Microbial	<100 CFU/g	No Growth

Approval

Quality Assurance – Print Name: <i>Terri Sorensen</i>	Quality Assurance – Sign Name: <i>TSorensen</i>	Date: <i>04/19/14</i>
--	--	--------------------------

Bovine albumin is derived from raw bovine serum processed from blood collected from healthy animals of U.S. Origin.