

## Sterile Purified Water, USP Grade

Product Number:	WPW-USP-1X6	Lot Number:	20712161
		Date of Manufacture:	12/2016
		Expiration Date:	12/2019
		Storage Temperature:	2 - 30° C

**BASIS FOR RELEASE**

CHEMICAL ANALYSIS			
Test	Method	Specification	Result
Appearance	Visual	Clear Liquid	Clear Liquid
Conductivity	USP <645>	<5.0 µS/cm	0.61 µS/cm
Sterility	USP <71>	No Growth	No Growth
TOC (before packaging)	USP <643>	<500 ppb	10.7 ppb

### STATEMENT FOR USE

This product is for further manufacturing, laboratory, or research use only.  
Not intended for direct human or veterinary parenteral administration or where prohibited by law.

**Approval**

Quality Assurance- Print Name:	Quality Assurance- Sign Name:	Date:
<b>Michelle Perea</b>	<i>Michelle Perea</i>	27 Dec 2016