

Sodium Chloride 0.9% (Normal Saline), USP Grade

Product Number: SAL-090-1X6

Lot Number: 21502171

Date of Manufacture: 02/2017

Expiration Date: 08/2018

Storage Temperature: 15 - 30 ° C


BASIS FOR RELEASE

CHEMICAL ANALYSIS			
Test	Method	Specification	Result
Appearance	Visual	Clear Liquid	Clear Liquid
Endotoxin	USP <85>	<0.25 EU/mL	< 0.005 EU/mL
Particulate Matter	USP <788>	≤ 3 particles per mL that are ≥ 25µm in size	0 particles/mL
Particulate Matter	USP <788>	≤ 25 particles per mL that are ≥ 10µm in size	0 particles/mL
pH	USP <791>	4.5 - 7.0	5.73
Sterility	USP <71>	No Growth	No Growth

STATEMENT FOR USE

This product is for further manufacturing, laboratory, or research use only.
Not intended for direct human or veterinary parenteral administration or where prohibited by law.

Approval

Quality Assurance- Print Name: Michelle Perea	Quality Assurance- Sign Name: 	Date: 07 MAR 2017
---------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------	-----------------------------