

## Non-Sterile Purified Water, USP Grade

Product Number: WPW-USN-2XL

Lot Number: 30401171

Date of Manufacture: 01/2017

Expiration Date: 01/2020

Storage Temperature: 2 - 30° C

**BASIS FOR RELEASE**

CHEMICAL ANALYSIS			
Test	Method	Specification	Result
Appearance	Visual	Clear Liquid	Clear Liquid
Conductivity	USP <645>	<5.0 µS/cm	0.63 µS/cm
TOC (before packaging)	USP <643>	<500 ppb	35.45 ppb

### STATEMENT FOR USE

This product is for further manufacturing, laboratory, or research use only.  
Not intended for direct human or veterinary parenteral administration or where prohibited by law.

#### Approval

Quality Assurance- Print Name: <b>Susan R. Riggs</b>	Quality Assurance- Sign Name: <i>Susan R. Riggs</i>	Date: <i>22 Feb 2017</i>
---	--	-----------------------------