

## Non-Sterile Purified Water, USP Grade

Product Number: WPW-USN-20L

Lot Number: 31506171

Date of Manufacture: 06/2017

Expiration Date: 06/2020

Storage Temperature: 2 - 30° C

**BASIS FOR RELEASE**

CHEMICAL ANALYSIS			
Test	Method	Specification	Result
Appearance	Visual	Clear Liquid	Clear Liquid
Conductivity	USP <645>	<5.0 µS/cm	0.69 µS/cm
TOC (before packaging)	USP <643>	<500 ppb	24.16 ppb

### STATEMENT FOR USE

This product is for further manufacturing, laboratory, or research use only.  
Not intended for direct human or veterinary parenteral administration or where prohibited by law.

#### Approval

Quality Assurance- Print Name:	Quality Assurance- Sign Name:	Date:
<b>Michelle Perea</b>	<i>Michelle Perea</i>	16 Jun 2017