

# Certificate of Analysis

## Purified Water, USP Sterile Grade

Product Number: WPW-USP-2XL

Lot Number: 31611182

Date of Manufacture: 11/2018

Expiration Date: 11/2021

Storage Temperature: 2 - 30° C

### BASIS FOR RELEASE

#### CHEMICAL ANALYSIS

Test	Method	Specification	Result
Appearance	Visual	Clear and Colorless Liquid	Clear and Colorless Liquid
Conductivity (Packaged)	USP <645>	≤ 5.0 μS/cm	0.42 μS/cm
Sterility	USP <71>	No Growth	No Growth
TOC (Before Packaging)	USP <643>	< 500 ppb	14.05 ppb

This lot has been produced using Current Good Manufacturing Practices (cGMP). The solution meets all requirements for residual solvents USP <467>.


This product is not of animal origin and does not contain, make use of or involve animal components during its manufacturing process.

No antimicrobial or other substance has been added.

### STATEMENT FOR USE

This product is for further manufacturing, laboratory, or research use only.  
 Not intended for direct human or veterinary parenteral administration or where prohibited by law.

#### Approval

Quality Assurance- Print Name: <b>Charles Burt</b>	Quality Assurance- Sign Name: 	Date: 06 Dec 2018
---	--	----------------------