

## Non-Sterile Purified Water, USP Grade

Product Number: WPW-USN-2XL	Lot Number: 33107171
	Date of Manufacture: 07/2017
	Expiration Date: 07/2020
	Storage Temperature: 2 - 30° C


**BASIS FOR RELEASE**

CHEMICAL ANALYSIS			
Test	Method	Specification	Result
Appearance	Visual	Clear Liquid	Clear Liquid
Conductivity	USP <645>	<5.0 µS/cm	0.55 µS/cm
TOC (before packaging)	USP <643>	<500 ppb	27.85 ppb

### STATEMENT FOR USE

This product is for further manufacturing, laboratory, or research use only.  
Not intended for direct human or veterinary parenteral administration or where prohibited by law.

**Approval**

Quality Assurance- Print Name: <b>Emily Sorensen</b>	Quality Assurance- Sign Name: 	Date: <b>01 Aug 2017</b>
---	--	-----------------------------