

Lipogro

Cholesterol – Lipoprotein Complex, 0.1 micron filtered; 100mL

Product number: RMB-001-201 Lot Number: 41115BCHD-R

Date of Manufacture: 4/2011 Expiration date: 4/2013

Storage temperature: 2 – 8°C

Not for human use

Basis for Release

Test / Method	Specification	Result
Cholesterol Concentration	9.00 – 11.00 g/L	9.78 g/L
Total Protein	13.0-18.0 g/L	14.0 g/L
pH	7.0 – 8.4	7.6
Endotoxin	< 50 EU/mL	35.7 EU/mL
Microbial Content	None Detected	Pass
Mycoplasma	None Detected	Pass
Virus	None Detected	Pass
Appearance	Clear yellow solution, visibly free of particulates	Pass
Filtration	0.1 µm	Pass

Approval

Quality Assurance – Print Name: <i>William D. Pfeffer</i>	Quality Assurance – Sign Name: <i>William D. Pfeffer</i>	Date: <i>12/14/2012</i>
--	---	----------------------------

Lipogro is derived from bovine serum, processed from blood collected from healthy animals. Characterization assays are performed on various lots of this product, and will be provided upon request.